



PRODUCT SHEET


FACE PROTECTION SHIELD

REFERENCE:

Blayco PF-01

COMPOSITION:

PARTS	MATERIAL
Clear shield	400 µm PET
Grip tape	50% Polyester 50% Elastane
Protective foam	Polyurethane
Rivet	Acetal resin
Individual bag with self-closing	PE
Secondary packaging	Carton box



DESCRIPTION AND INDICATION OF USE:

Face Protection Shields are disposable protection devices for healthcare professionals, exposed to possible contamination by microorganisms spread by droplet. They serve to protect the facial area and associated mucous membranes (eyes, nose, mouth) against splashes of body fluids, aerosols and contaminating droplets.

Telic, S.A.U guarantees an adequate microbiological quality of its Blayco® Face Protection Shields PF-01, through good manufacturing practices certified as ISO13485 for the production of medical devices. If you want to disinfect the product or if its re-use is necessary, it is recommended to clean it with disinfectants such as 62-71% ethanol, 0.5% hydrogen peroxide or 0.1% sodium hypochlorite.

Product classified as Personal Protective Equipment Cat. 1. Non-sterile.

QUANTITY PER BOX: 50 units

STORAGE CONDITIONS

Keep in a cool and dry place, at a temperature of 5 ° to 40 ° C.